



Simtrak™ Mobility

A decade-long background in professional and amateur sport medical data collection, analysis, and reporting.

I N S T A L L the Simtrak™ Help Guide series

Simtrak™ Mobility 9 Installation Instructions

Thank you for purchasing Simtrak™ Mobility 9. All of the necessary components to get started are on the CD-ROM.

Folders included on the CD-ROM

1. Documentation – This document, version 9.0 release notes, setup instructions, license agreement, and other technical documents.
2. Remove8 – REMOVE80.EXE. Program used to clean out release 8.0 files before update is run for a new release.
3. Installation – Installation for server or PC. **New installations only. Do not use this installation for updates (new releases).**
4. Upgrade – Installation for updating from version 8.0. **Use this installation for updates (new releases).**
5. Workstation – Installation for a workstation connected to a server where Simtrak™ has already been installed. Will only install run-time DLL's and create an empty simtrak9 folder on the workstation. Shortcuts to the server will then have to be created manually.

The instructions on the other side of this page provide the steps for installing the software to a new environment, where there has not been a Simtrak™ Mobility installation previously.

If there is a previous installation of Simtrak™ Mobility at your location, then put aside this document and use the Help Guide series document entitled UPDATE.

Installing The Software, for new users (first time installation), on single computer or on network server

1. Put the product CD-ROM in the CD-ROM drive (for network install use a workstation, not a server).
2. Click on the Windows **Start** button and choose **Run**.
3. Click on the **Browse** button.
4. Select the CD-ROM drive from the **Look In** pull-down.
5. Double-click on the folder named "Installation".
6. Click on the icon labeled **SimtrakMobility9-Installation-SS** (or **-SM, -PM, -PS,**) and then on the **Open** button.
7. Click on the **OK** button.
8. A message will be displayed that the installation is preparing for install.
If a message comes up that the Windows Installer is an older version, click on **OK** to continue. If your version of Windows Installer is not current enough, the installation will install a current version and notify you that you must re-boot your computer to continue. Select OK and let the computer re-boot. The installation process will then continue.
9. The **Welcome** screen will be displayed, with copyright information. Click the **Next** button.
10. The **Software License Agreement** screen will display. Accept the terms. Click the **Next** button.
11. The **Customer Information** screen will display. Verify the information is correct. Click the **Next** button.
12. The **Destination Folder** screen will display. The default destination is **c:\simtrak9**, the location generally used for installations on a single computer. To change the destination folder if you like, click on the **Change** button and select a different folder (for a network install, you must select the designated server location ... mapped drive). With the correct folder selected, click the **Next** button.
13. The **Ready To Install** screen will be displayed. Click the **Install** button.
14. The software will be installed.
If an error message is displayed saying "1907 Could not register font", click Ignore (this is not a problem to consider), and the installation will continue.
15. When the installation is finished, the **Completed** screen will display. A **Simtrak™ Mobility 9** icon will be on the desktop that was used for the installation. Click the **Finish** button. You are done.

For Network Installations – setting up workstations (after software is installed on server)

NOTE: All users need administrative rights to the simtrak root folder and its sub folders (WRITE and DELETE admin permissions in addition to READ and EXECUTE user permissions). Otherwise, will get "Resource file not valid ..." and "Cannot update cursor ..." error messages.

1. Place the Simtrak™ Mobility CD in the CD-ROM drives of each of the workstations that were not used for the Simtrak™ Mobility software installation.
2. At each of these workstations, double-click on the folder named "Workstation" on the CD.
3. Run the file named **SimtrakMobility91 – Workstation**. This step will place certain DLL's and an empty folder named simtrak9 on the workstation C drive. (*The workstation used for the initial server installation will not have this empty folder created on that workstation hard drive.*) You can remove this empty simtrak9 folder from the workstation C drive; it is not needed.
4. Remove the CD from the drive.
5. Create a shortcut on the workstation desktop (right click, New/Shortcut). Browse to the simtrak9 folder (or other folder you have designated) on the server, select simtrak9.exe, click OK, then Next. Enter the name for the icon, and click Finish. If you want to verify the settings, right click on the icon and choose Properties. The Target should be the simtrak9 (or other) folder and the simtrak9.exe file; the Start folder should be the simtrak9 (or other) folder.
6. Do this for each workstation that will be using Simtrak™ Mobility on the network.